

TDRKH specific antibody



Catalog Number: 115930

Product name

TDRKH specific antibody

IP: 1:2000-1:20000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

TDRKH specific Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in Hela cells. Positive WB detected in mouse testis tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 67kd

Validations




Preparation

This antibody was obtained by immunization of TDRKH specific recombinant protein (Accession Number: BC032690). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115930(TDRKH antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115930(TDRKH antibody) at dilution of 1:100 (under 40x lens)




TDRKH specific antibody

Catalog Number: 115930

Immunofluorescent analysis of HeLa cells, using TDRKH antibody Catalog No:115930 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IP Result of anti-TDRKH (IP:Catalog No:115930, 4ug; Detection:Catalog No:115930 1:7000) with mouse brain tissue lysate 4000ug.



IgG (red).