

# Anti-PTN Rabbit antibody



Catalog Number: 166930

## Product name

Anti-PTN Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PTN Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human PTN

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000-1/5000

IP: 1/20

## Validations



Western blot detection of PTN in C6,Rat Brain cell lysates using PTN Rabbit pAb(1:1000 diluted).Predicted band size:19kDa.Observed band size:19kDa.