

# Vorsetuzumab reference antibody



Catalog Number: 176930

## Product name

Vorsetuzumab reference antibody

## Specificity

Human TNFSF7 / CD27L / CD70

## Antibody description

Anti-TNFSF7 / CD27L / CD70 Antibody  
(Vorsetuzumab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-TNFSF7 / CD27L / CD70 Antibody (vorsetuzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-TNFSF7 / CD27L / CD70 Antibody (vorsetuzumab) is 99.29%, determined by SEC-HPLC.



## ELISA

Immobilized humanCD70-His at 2 µg/mL (30 µL/well) can bind Vorsetuzumab(P98623) with the EC50 of 0.001164 µg/mL.