

Human VCP (His tag) recombinant protein



Catalog Number: 518930

General Information

Gene Name Synonym

Transitional Endoplasmic Reticulum ATPase, TER ATPase, 15S Mg(2+)-ATPase p97 Subunit, Valosin-Containing Protein, VCP

Protein Construction

A DNA sequence encoding the Human VCP Ala2-Asn589 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM TrisHCl pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human VCP (His tag) recombinant protein