

GTF2F1 antibody

Catalog Number: 111230

Product name

GTF2F1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GTF2F1 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of GTF2F1 recombinant protein (Accession Number: NM_002096). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111230(GTF2F1 antibody) at dilution of 1:1000



IP Result of anti-GTF2F1 (IP:Catalog No:111230, 4ug; Detection:Catalog No:111230 1:300) with K-562 cells lysate 4000ug.