

Anti-PREP antibody



Catalog Number: 142230

Product name

Anti-PREP antibody

Immunogen

[Mouse Prolyl endopeptidase / PREP \(His Tag\) recombinant protein](#)

Specificity

Mouse Prolyl endopeptidase / PREP

Antibody description

Rabbit Polyclonal to Mouse PREP

Preparation

Produced in rabbits immunized with purified, recombinant Mouse Prolyl endopeptidase / PREP . Prolyl endopeptidase / PREP specific IgG was purified by Mouse Prolyl endopeptidase / PREP affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 1-4 μ l/mg of lysate

Validations



Lanes	A	B
Sample (whole cell lysate)	MCF-7	293T
Sample quantity	0.5 mg	
IP antibody quantity	4 μ L	
Protein G agarose	15 μ L of 50% Protein G Agarose	
Gel	13% SDS-PAGE reducing gel	
Primary antibody	His- α PREP antibody at 10 μ g/ml	
Secondary antibody	Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:2000 dilution	

Prolyl endopeptidase / PREP Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation