

# FYB Monoclonal Antibody



Catalog Number: 166230

## Product name

FYB Monoclonal Antibody

## Specificity

Human

## Antibody description

Monoclonal antibody to FYB

## Preparation

Antigen: Recombinant protein of human FYB

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

WB 1:200 - 1:1000

## Validations



Western blot - FYB Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using FYB antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.