

# TRAF3 antibody

Catalog Number: 116230

## Product name

TRAF3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TRAF3 Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in HEK-293 cells, COLO 320 cells, HeLa cells, Raji cells. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of TRAF3 recombinant protein (Accession Number: XM\_017021618). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


## Validations



WB results of Catalog No:(anti TRAF3).



IP Result of anti-TRAF3 (IP:Catalog No:116230, 4ug; Detection:Catalog No:116230 1:800) with



Raji cells lysate 2000ug.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116230(TRAF3 Antibody) at dilution of 1:600