

ASIC4 antibody



Catalog Number: 108230

Product name

ASIC4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ASIC4 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Observed molecular weight by Western-blot: 70kd

Preparation

This antibody was obtained by immunization of ASIC4 recombinant protein (Accession Number: BC010439). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:108230(ACCN4 antibody) at dilution of 1:400