

# COL17A1 Monoclonal Antibody



Catalog Number: 160630

## Product name

COL17A1 Monoclonal Antibody

IHC 1:50 - 1:200

IF 1:50 - 1:200

IP 1:20 - 1:50

## Specificity

Human, Mouse, Rat

## Antibody description

Monoclonal antibody to COL17A1

## Preparation

Antigen: Recombinant protein of human COL17A1

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC, IF, IP

## Dilutions


WB 1:500 - 1:2000

## Validations



Immunohistochemistry - COL17A1 Monoclonal Antibody

Immunohistochemistry of paraffin-embedded mouse skin using COL17A1 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - COL17A1 Monoclonal Antibody

Immunofluorescence analysis of HUVEC cells using COL17A1 antibody .