

# TIRAP Polyclonal Antibody



Catalog Number: 163630

## Product name

TIRAP Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to TIRAP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human TIRAP (NP\_001034750.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - TIRAP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using TIRAP antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - TIRAP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using TIRAP antibody at dilution of 1:200 (40x lens).