

# TNFAIP1 antibody



Catalog Number: 107630

## Product name

TNFAIP1 antibody

## Specificity

Human; other species not tested.

## Antibody description

TNFAIP1 Mouse Monoclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of TNFAIP1 recombinant protein (Accession Number: NM\_021137). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgM

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107630(TNFAIP1 Antibody) at dilution of 1:500