

Anti-TAS2R8 antibody



Catalog Number: 176330

Product name

Anti-TAS2R8 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TAS2R8

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tas2r8. aa range:45-94

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.: Western Blot analysis of various cells using T2R8 Polyclonal Antibody