

# FMO5 antibody

Catalog Number: 110730

## Product name

FMO5 antibody

## Immunogen

[Human FMO5 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FMO5 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse lung tissue. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of FMO5 recombinant protein (Accession Number: NM\_001144830). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110730(FMO5 antibody) at dilution of 1:1200



IP Result of anti-FMO5 (IP:Catalog No:110730, 4ug; Detection:Catalog No:110730 1:500) with mouse lung tissue lysate 4000ug.