

RNASEH2B antibody



Catalog Number: 114730

Product name

RNASEH2B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RNASEH2B Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HeLa cells, MCF7 cells, mouse lung tissue. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of RNASEH2B recombinant protein (Accession Number: BC036744). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114730(RNASEH2B antibody) at dilution of 1:500