

DXd Monoclonal Antibody



Catalog Number: 176730

Product name

DXd Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

DXd conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4), lyophilized. **Reconstitution: Adding sterile water, prepare a stock solution of 2.0 mg/ml.**

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1, κ


Applications

ELISA

Dilutions


ELISA: 1:5,000-1:20,000

Validations



ELISA Binding Assay of Anti-DXd Antibody to different proteins.

ELISA analysis of Anti-DXd Monoclonal Antibody with Dxd-BSA, BSA and Dxd ADC. Anti-DXd Monoclonal Antibody can bind to Dxd-BSA or Dxd ADC but not carrier protein.



ELISA Binding Assay of Anti-DXd Antibody

Immobilized Dxd ADC at 2 μ g/mL (100 μ L/well) can bind anti-DXd mouse monoclonal antibody.