

# VDAC2 antibody

Catalog Number: 116730

## Product name

VDAC2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VDAC2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human prostate tissue. Positive WB detected in human heart tissue, human liver tissue. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 32kd

## Preparation

This antibody was obtained by immunization of VDAC2 recombinant protein (Accession Number: NM\_001184823). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116730(VDAC2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:116730(VDAC2 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-VDAC2 (IP:Catalog No:116730,

# VDAC2 antibody



Catalog Number: 116730

---

4ug; Detection:Catalog No:116730 1:500) with mouse heart tissue lysate 3200ug.