

CASP12 Polyclonal Antibody



Catalog Number: 160040

Product name

CASP12 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CASP12

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CASP12 (NP_001177945.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - CASP12 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CASP12 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - CASP12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using CASP12 antibody at dilution



of 1:100 (20x lens).

CASP12 Polyclonal Antibody



Catalog Number: 160040

Immunohistochemistry - CASP12 Polyclonal
Antibody

Immunohistochemistry of paraffin-embedded
human stomach using CASP12 antibody at dilution
of 1:100 (40x lens).