

HOXD13 antibody



Catalog Number: 111440

Product name

HOXD13 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

HOXD13 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human hepat°Cirrhosis tissue, human liver cancer tissue. Positive WB detected in A549 cells, HEK-293 cells, HepG2 cells, Jurkat cells, mouse brain tissue, mouse lung tissue, mouse small intestine tissue, rat brain tissue, Recombinant protein. Positive IP detected in A549 cells. Observed molecular weight by Western-blot: 36kd

Validations




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000523). Purification method: Antigen affinity purified.

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111440(HOXD13 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



ELISA, WB, IHC, IP

Immunohistochemical of paraffin-embedded human brain using Catalog No:111440(HOXD13 antibody) at dilution of 1:100 (under 10x lens)

Dilutions

HOXD13 antibody

Catalog Number: 111440



Immunohistochemical of paraffin-embedded human brain using Catalog No:111440(HOXD13 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-HOXD13 (IP:Catalog No:111440, 4ug; Detection:Catalog No:111440 1:500) with A549 cells lysate 1600ug.