

Anti-HLA-DQB1 antibody

Catalog Number: 176440

Product name

Anti-HLA-DQB1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to HLA-DQB1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human hla-dqb1/hla-dqb2. aa range:131-180

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations





Fig1.; Western Blot analysis of PC12 cells using HLA-DQB1/2 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at



1:20000

Fig2.; Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:100