

Anti-Brk/PTK6 (7D6) Mouse antibody



Catalog Number: 168440

Product name

Anti-Brk/PTK6 (7D6) Mouse antibody

Specificity

Human

Antibody description

Brk/PTK6 (7D6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PTK6 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, IP

Dilutions


WB: 1/500 - 1/2000

ELISA: 1/10000

Validations




Western blot analysis using PTK6 mouse mAb against HeLa (1), A549 (2) and MCF-7 (3) cell



lysate.

Immunoprecipitation analysis of HeLa cell lysates



using Brk mouse mAb.

Western blot detection of Brk in K562, THP-1, MCF7, A549 and HeLa cell lysates using Brk mouse mAb (1:1000 diluted). Predicted band size: 52KDa. Observed band size: 52KDa.