

# Mouse Contactin 3/CNTN3 (His Tag) recombinant protein



Catalog Number: 500840

## General Information

### Gene Name Synonym

Brain-derived immunoglobulin superfamily protein 1; Plasmacytoma-associated neuronal glycoprotein

### Protein Construction

A DNA sequence encoding the mouse CNTN3 (NP\_032805.2) (Met1-Gly 1001) was expressed, with a C-terminal polyhistidine tag.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu 20

### Molecular Mass

The recombinant mouse CNTN3 consists of 993

amino acids and predicts a molecular mass of 110 KDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmCNTN3 is approximately 140-150 KDa as a result of glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse Contactin 3 / CNTN3 Protein (His Tag) SDS-PAGE