

Anti-phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) Rabbit antibody



Catalog Number: 167840

Product name

Anti-phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

phospho-Ezrin (Thr567)/ Radixin (Thr564)/ Moesin (Thr558) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr567 of human Ezrin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of Ezrin (phospho-Thr567)/ Radixin (phospho-Thr564)/ Moesin (phospho-Thr558) in HeLa, CHO-K1, C6 using Ezrin (phospho-Thr567)/ Radixin (phospho-Thr564)/ Moesin (phospho-Thr558) Rabbit mAb(1:1000 diluted)