

# DYNLL2 antibody



Catalog Number: 110140

---

## Product name

DYNLL2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DYNLL2 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, A431 cells, human brain tissue. Observed molecular weight by Western-blot: 10kd

## Preparation

This antibody was obtained by immunization of DYNLL2 recombinant protein (Accession Number: NM\_080677). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110140(DYNLL2 antibody) at dilution of 1:400