

NFAM1 antibody

Catalog Number: 113140

Product name

NFAM1 antibody

Specificity

Human; other species not tested.

Antibody description

NFAM1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of NFAM1 recombinant protein (Accession Number: NM_145912). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113140(NFAM1 Antibody) at dilution of 1:350



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:113140(NFAM1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:113140(NFAM1 Antibody) at dilution of 1:50 (under 40x lens)