

# NOSTRIN Polyclonal Antibody



Catalog Number: 164140

## Product name

NOSTRIN Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to NOSTRIN

## Preparation

Antigen: A synthetic peptide of human NOSTRIN

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - NOSTRIN Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using NOSTRIN antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.