

# SF3B5 antibody



Catalog Number: 115140

## Product name

SF3B5 antibody

## Immunogen

[Human SF3B5 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SF3B5 Rabbit Polyclonal antibody. Positive WB detected in human adrenal gland tissue, HeLa cells, mouse eye tissue, SH-SY5Y cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 10-12 kDa

## Preparation

This antibody was obtained by immunization of SF3B5 recombinant protein (Accession Number: NM\_031287). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:115140(SF3B5 antibody) at dilution of 1:400



Immunofluorescent analysis of HeLa cells, using SF3B5 antibody Catalog No:115140 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).