

FBX15 antibody



Catalog Number: 110540

Product name

FBX15 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FBX15 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, mouse ovary tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of FBX15 recombinant protein (Accession Number: XM_017025604). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:110540(FBXO15 antibody) at dilution of 1:1500