

MBTD1 antibody



Catalog Number: 112540

Product name

MBTD1 antibody

Specificity

Human; other species not tested.

Antibody description

MBTD1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HL-60 cells, K-562 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 65-71 kDa

Preparation

This antibody was obtained by immunization of MBTD1 recombinant protein (Accession Number: BC101736). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112540(MBTD1 Antibody) at dilution of 1:600