

# SPANXB1 antibody



Catalog Number: 115540

## Product name

SPANXB1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SPANXB1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of SPANXB1 recombinant protein (Accession Number: NM\_145664). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115540(SPANXB1 Antibody) at dilution of 1:300