

# Anti-GC1q R Rabbit antibody



Catalog Number: 167540

## Product name

Anti-GC1q R Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GC1q R Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human C1QBP

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000-1/5000

IP: 1/20

## Validations



Western blot detection of GC1q R in K562,Rat Brain,C6,3T3,Hela cell lysates using GC1q R Rabbit pAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:31kDa.