

KHDRBS2 Polyclonal Antibody



Catalog Number: 162940

Product name

KHDRBS2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to KHDRBS2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human KHDRBS2 (NP_689901.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - KHDRBS2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using KHDRBS2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.