

Anti-CYP3A7-CYP3A51P antibody



Catalog Number: 175940

Product name

Anti-CYP3A7-CYP3A51P antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CYP3A7-CYP3A51P

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cytochrome p450 3a7. aa range:231-280

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293T cells using CYP3A7 Polyclonal Antibody diluted at 1: 500