

TEKT5 antibody



Catalog Number: 115940

Product name

TEKT5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TEKT5 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, HepG2 cells, mouse brain tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of TEKT5 recombinant protein (Accession Number: NM_144674). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115940(TEKT5 Antibody) at dilution of 1:1000