

Anti-CD79a Rabbit antibody



Catalog Number: 168940

Product name

Anti-CD79a Rabbit antibody

Specificity

Human

Antibody description

CD79a Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human CD79a

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/2000

IHC: 1/20-1/50

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of CD79a in Jurkat cell lysates using CD79a Rabbit mAb(1:500 diluted).Predicted band size:25kDa.Observed band size:45-55kDa.