

Human PRMT6 (His & FLAG Tag) recombinant protein



Catalog Number: 501240

General Information

Gene Name Synonym

Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6; Histone-arginine N-methyltransferase PRMT6

Protein Construction

A DNA sequence encoding the human PRMT6 (NP_060607.2) (Met 1-Asp 375) was expressed, fused with a polyhistidine tag at the C-terminus and the flag tag at the N-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human PRMT6 consists of 394 amino acids and predicts a molecular mass of 44.4 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh PRMT6 is approximately 43-46 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human PRMT6 / HRMT1L6 Protein (His & FLAG Tag) SDS-PAGE