

Mouse EpCAM/TROP-1/TACSTD1 (His Tag) recombinant protein



Catalog Number: 501640

General Information

Gene Name Synonym

Epithelial glycoprotein 314; Protein 289A; Tumor-associated calcium signal transducer 1

Protein Construction

A DNA sequence encoding the extracellular domain of mouse EpCAM (NP_032558.2) (Met 1-Thr 266) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 24

Molecular Mass

The recombinant mouse EpCAM consists of 254

amino acids and has a predicted molecular mass of 29 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rmEpCAM is approximately 37 kDa due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse EpCAM / TROP-1 / TACSTD1 Protein (His Tag) SDS-PAGE