

Anti-TYMS antibody



Catalog Number: 142640

Product name

Anti-TYMS antibody

Specificity

Human TYMS

Antibody description

Rabbit Polyclonal to Human TYMS

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human TYMS, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions


WB: 5-20 µg/ml

IP: 4-8 µl/mg of lysate

Validations



Human TYMS Western blot (WB) 14654



Human TYMS Immunoprecipitation(IP) 14655