

CHIC2 Polyclonal Antibody



Catalog Number: 165640

Product name

CHIC2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CHIC2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human CHIC2 (NP_036242.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - CHIC2 Polyclonal Antibody

Western blot analysis of extracts of mouse testis, using CHIC2 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - CHIC2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CHIC2 antibody at dilution



of 1:100 (40x lens).


CHIC2 Polyclonal Antibody



Catalog Number: 165640

Immunohistochemistry - CHIC2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using CHIC2 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CHIC2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using CHIC2 antibody at dilution of 1:100 (40x lens).