

# Anti-LTA antibody



Catalog Number: 176640

## Product name

Anti-LTA antibody

## Ig Type

Mouse IgG2b

## Specificity

Human

## Applications

WB

## Antibody description

Mouse monoclonal antibody to LTA

## Dilutions

WB: 1:2,000

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Validations




Fig1: Western blot analysis on recombinant LT-alpha lysates using anti-LT alpha monoclonal antibody.