

Anti-LTA antibody

Catalog Number: 176640

Product name

Anti-LTA antibody

Specificity

Human

Antibody description

Mouse monoclonal antibody to LTA

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

WB

Dilutions

WB: 1:2,000

Validations




Fig1: Western blot analysis on recombinant LT-alpha lysates using anti-LT alpha monoclonal antibody.