

# PURB antibody



Catalog Number: 114340

## Product name

PURB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PURB Rabbit Polyclonal antibody. Positive WB detected in C2C12 cell, mouse liver tissue. Observed molecular weight by Western-blot: 33-40 kDa

## Preparation

This antibody was obtained by immunization of PURB recombinant protein (Accession Number: NM\_033224). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



C2C12 cell were subjected to SDS PAGE followed by western blot with Catalog No:114340(PURB Antibody) at dilution of 1:600