

Anti-Cardiac Troponin I (1D7) Mouse antibody



Catalog Number: 168340

Product name

Anti-Cardiac Troponin I (1D7) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

Cardiac Troponin I (1D7) Mouse monoclonal antibody

Preparation

Antigen: Synthetic peptides corresponding to the sequence of human Cardiac Troponin I.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, ELISA

Dilutions


ELISA: 1/10000

;WB: 1/1000

Validations



sELISA standard Curve for Cardiac Troponin I: Capture Antibody Mouse mAb to Cardiac Troponin I at 2ug/ml and 168340(1D7-D10) was used for



detecting.

Western blot analysis of extracts from Rat and mouse heart tissue using Troponin I (168340,dilution 1:1000) Mouse mAb.Predicted band size:28kDa.Observed band size:28kDa.