

Claudin 23 antibody



Catalog Number: 109340

Product name

Claudin 23 antibody

Specificity

Human; other species not tested.

Antibody description

Claudin 23 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of Claudin 23 recombinant protein (Accession Number: NM_194284). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:109340(CLDN23 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:109340(CLDN23 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:109340(CLDN23 Antibody) at dilution of 1:50 (under 40x lens)