

# FNTA antibody



Catalog Number: 110740

## Product name

FNTA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FNTA Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, 3T3-L1 cells, A431 cells, Jurkat cells. Observed molecular weight by Western-blot: 44 kDa

## Preparation

This antibody was obtained by immunization of FNTA recombinant protein (Accession Number: NM\_002027). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110740(FNTA antibody) at dilution of 1:500