

Anti-HA antibody

Catalog Number: 141740

Product name

Anti-HA antibody

Immunogen

[Influenza B HA \(His Tag\) recombinant protein](#)

Specificity

Influenza B virus (B/Florida/4/2006)
Hemagglutinin / HA

Antibody description

Rabbit Polyclonal to Influenza B HA

Preparation

Produced in rabbits immunized with purified, recombinant Influenza B virus (B/Florida/4/2006) Hemagglutinin / HA . Influenza B virus (B/Florida/4/2006) Hemagglutinin / HA specific IgG was purified by Influenza B virus (B/Florida/4/2006) Hemagglutinin / HA affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ELISA, IHC-P, FCM, ICC/IF, IF, IP


Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza B virus (B/Florida/4/2006) Hemagglutinin / HA.

Validations



Items	Lanes	A	B
Sample (Recombinant protein)	B-F-HA (Cat#11053-V08H)		
Sample Volume (ng/lane)	50	10	
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-5 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			

Influenza B HA Western blot (WB) 14705