

STT3B antibody



Catalog Number: 115740

Product name

STT3B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STT3B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 94 kDa

Preparation

This antibody was obtained by immunization of STT3B recombinant protein (Accession Number: BC015880). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115740(STT3B antibody) at dilution of 1:2000