

SAMD11 antibody



Catalog Number: 115050

Product name

SAMD11 antibody

Specificity

Human; other species not tested.

Antibody description

SAMD11 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of SAMD11 recombinant protein (Accession Number: BC033213). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115050(SAMD11 Antibody) at dilution of 1:600