

Recombinant AMOL1 protein (His tag)



Catalog Number: 538450

Gene Name Synonym

Class B basic helix-loop-helix protein 25; bHLHb25; Helix-loop-helix protein HEIR-1; ID-like protein inhibitor HLH 1R21; Inhibitor of DNA binding 3; Inhibitor of differentiation 3

Protein Construction

DNA encoding 432-777 aa (Q8IY63) were fused with C-terminal 6His tag.

Organism

Human

Expression Host

E. coli

Activity

Not tested.

Purity

>85%, by SDS-PAGE with Coomassie Brilliant Blue staining.

Endotoxin

Please contact the lab for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Molecular Mass

Predictes a molecular mass of 41.90 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 50 kDa band.

Formulation

Lyophilized from sterile 50 mM HEPES, 250 mM NaCl, pH 7.4


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

SDS-PAGE



Recombinant AMOL1 protein (His tag) SDS-PAGE

M: Marker; Lane 1: AMOL1; Lane 2: BSA