

# WBSCR16 antibody



Catalog Number: 116850

## Product name

WBSCR16 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

WBSCR16 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of WBSCR16 recombinant protein (Accession Number: NM\_148842). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116850(WBSCR16 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:116850(WBSCR16 antibody) at dilution of 1:50



(under 40x lens)

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116850(WBSCR16 Antibody) at dilution of

# WBSCR16 antibody



Catalog Number: 116850

---

1:600